



Liverpool Surgical Supplies

Unit 1/30 Heathcote Rd, Moorebank NSW 2170

Ph: (02) 9602 0069 Fax: (02) 9601 7870

Email: karync@liverpoolsurgicalsupplies.com.au

Website: [www.hospitalathome.com.au/
product-category/compression-stockings/](http://www.hospitalathome.com.au/product-category/compression-stockings/)

Consult is via appointment only.

Please call to book your appointment
on (02) 9602 0069



Liverpool Surgical Supplies is a division of
PEGASUS
HEALTH GROUP



Specialising in
ready to wear
surgical compression
garments



DVA and NDIS approved

Garments are recommended for a variety of reasons

- Tired aching limbs
- Swollen legs or ankles
- Poor circulation
- Travel
- Spider veins
- Pregnancy
- DVT (deep vein thrombosis)
- Leg ulcers
- Arm oedema

Graduated Compression Stockings, Sleeves and Wraps

Class 1- Minor foot, ankle and leg swelling, tired aching legs, pregnancy, travel

Class 2 - Swollen legs or arms, moderate varicose veins, post vein surgery/sclerotherapy, superficial thrombophlebitis,

Class 3 - Pronounced varicose veins, leg ulcers, constant swelling in ankles/legs, deep vein thrombosis, lymphatic oedema

Farrow Wrap - Lymphoedema and wound management

Preventative Measures

- Regular exercise
- Lose excess weight
- Avoid standing or sitting for long periods of time
- Avoid extreme heat i.e. hot baths, overexposure to sun
- Avoid lifting or carrying heavy loads
- Avoid constrictive clothing or footwear
- Always sustain a diet high in fibre
- Always drink sufficient fluids, preferably 2 litres of water per day

Whilst compression garments may be worn for a variety of conditions, we highly recommend that appropriate medical advice has been sought prior to purchase of your compression garment.



Order For:

Prescribing Dr: _____

Date: ____ / ____ / ____ Tel: _____


1. Select Brand

- JOBST VENOSAN MEDI
 IBICI JUZO SIGVARIS
 ANTI EMBOLISM
Anti-embolic stockings are intended for use during surgical procedures and on patients confined to bed. They are designed to prevent blood clots rather than active oedema management.

2. Select Compression

- Class 1 - 15-20mmHg
 Class 2 - 20-30mmHg
 Class 3 - 30-40mmHg
 Farrow Wrap - On Application


3. Indicate Style Required



4. Other Comments

Arm Circumferences		Leg Circumferences	
left	right	left	right
_____	H _____	_____	H _____
_____	G _____	_____	G _____
_____	F _____	_____	D _____
_____	D _____	_____	C _____
_____	C _____	_____	B _____

- H.** The length between G. and C. with arm slightly bent
G. Measure circumference 3cm from axilla
F. Measure circumference of mid upper-arm
D. Measure circumference of mid forearm
C. Measure circumference of wrist



- H.** Measure from the back of knee to floor or thigh to the floor depending on length required
G. Measure circumference of thigh at highest point
D. Measure circumference 1cm below knee
C. Measure circumference of calf at largest point
B. Measure circumference of ankle at smallest point

